

Measuring Dyskinesia in Parkinson's Clinical Trials

Thomas Mera, MS
Sr. Biomedical Engineer
February 21st, 2013
Webinar Series




Question 1:

1. How are PD motor complications evaluated in clinical trials?
2. What are the challenges with clinical trial dyskinesia endpoints, and how can they be improved?
3. How can home-based motion sensor dyskinesia assessment improve your clinical trials?

Motor Complications of Chronic Levodopa Therapy

1. Motor fluctuations
 - Alternate between therapy “off” and “on” states over dose cycles
2. Levodopa-induced dyskinesia (LID)
 - Involuntary, episodic, and irregular movements
 - Peak-dose most common

Advanced Stages



Keijsers, N. L., M. W. Horstink, et al. (2003). "Automatic assessment of levodopa-induced dyskinesias in daily life by neural networks." *Mov Disord* 18(1): 70-80.

Question 2:

1. How are PD motor complications evaluated in clinical trials?
2. What are the challenges with clinical trial dyskinesia endpoints, and how can they be improved?
3. How can home-based motion sensor dyskinesia assessment improve your clinical trials?

Clinical Trial Review

ClinicalTrials.gov

A service of the U.S. National Institutes of Health

Example: "Heart attack" AND "Los Angeles"

Search for studies:

Search

[Advanced Search](#) | [Help](#) | [Studies by Topic](#) | [Glossary](#)

[Find Studies](#) ▾

[About Clinical Studies](#) ▾

[Submit Studies](#) ▾

[Resources](#) ▾

[About This Site](#) ▾

[Home](#) > [Find Studies](#) > [Search Results](#)

[Text Size](#) ▾

377 studies found for: "parkinson's disease" AND "dyskinesia"

[Modify this search](#) | [How to Use Search Results](#)

List

By Topic

On a Map

Search Details

+ Show Display Options

 Download

 Subscribe to RSS

Include only open studies Exclude studies with unknown status

Rank	Status	Study
1	Not yet recruiting	Safety and Efficacy of AVP-923 in the Treatment of Levodopa-induced Dyskinesia in Parkinson's Disease Patients Conditions: Dyskinesia; Parkinson's Disease Interventions: Drug: AVP-923-45; Drug: Placebo
2	Recruiting	Open-label, Long-term Safety Extension Study of AFQ056 in Parkinson's Patients With L-dopa Induced Dyskinesias Conditions: Dyskinesias; Parkinson Disease; Movement Disorders; Parkinsonian Disorders; Anti-Dyskinesia Agents

Clinical Trial Endpoints

- Clinical Assessments
 - UPDRS, UDysRS, mAIMS, PDYS-26
 - Patient retrospective recall
- Patient diaries
 - Self assessment at home
 - 0.5-1 hr interval diary entries
- Body-worn motion sensors
 - Shift in research
 - Unconstrained continuous assessment at home
 - Transition to clinical trial use not trivial
 - Importance of quality assurance (e.g. FDA, ISO, CE, TGA)



Challenges with Clinical Trial Endpoints

- Resolution of clinical rating scales
 - Severity: 0-4 integer scoring
 - Temporal: snapshot of dyskinesia response
- Compliance of home diaries
 - Correlation between reported and actual compliance
 - Patient awareness of, understanding, and recognizing therapy states
- Costs
 - Clinician and patient time in clinic
 - Accuracy may affect statistical power

Motion Sensor LID Assessment: Clinical Validation Study

- Collaborators

- Michelle Burack, MD, PhD
- NIH-funded SBIR Phase I





UNIVERSITY of
ROCHESTER
MEDICAL CENTER

- Goals

1. Capture peak-dose dyskinesia over a levodopa dose using hand-worn motion sensors
2. Develop a scoring model to automatically rate dyskinesia
3. Determine whether a single motion sensor unit could accurately assess global dyskinesia

Methods: Study Preparation


- Off levodopa from previous night or end of dose
- A wireless motion sensor unit positioned on each hand
- Two discrete motor tasks:
 1. Arms resting
 2. Arms Extended
- Serial subtractions as distraction



Methods: Data Collection


- Two motor tasks at hours 0, 1, 2, and 3 after levodopa dose
- Motion sensor data were wirelessly streamed to a computer
- Video of task performance was recorded and later scored by two expert raters
 - modified-Abnormal Involuntary Movement Scale (m-AIMS)
 - 0 (none) to 4 (severe) global dyskinesia ratings
- Severity scoring models developed using sensor data and clinician global m-AIMS scores

Clinical Assessment





The time to reach peak-dose dyskinesia varied by subject

Symptom Feature Extraction






Dyskinesia Severity Scoring Model



Question 3:

1. How are PD motor complications evaluated in clinical trials?
2. What are the challenges with clinical trial dyskinesia endpoints, and how can they be improved?
3. How can home-based motion sensor dyskinesia assessment improve your clinical trials?


Clinical Trial Drug Comparison



Advantages For Your Clinical Trials

Motion sensor assessment during discrete tasks

- Clinical validation and quality assurance
- Home assessment kit
- Single motion sensor to assess global dyskinesia fluctuations
- Electronic formatting
 - Instant access, reports



Advantages For Your Clinical Trials

Motion sensor assessment during discrete tasks

- Clinical validation and quality assurance
- Home assessment kit
- Single motion sensor to assess global dyskinesia fluctuations
- Electronic formatting
 - Instant access, reports

Time	Rest Tremor	Postural Tremor	Finger Taps Speed	Finger Taps Amplitude	Finger Taps Rhythm	Dyskinesia
6:55 AM	3.9	3.4	2.6	2.5	2.3	0.0
6:57 AM	Drug Dose A					
7:28 AM	2.5	3.0	1.7	1.4	1.0	0.0
7:59 AM	0.5	1.9	1.8	1.5	1.2	1.3
8:30 AM	0.3	0.9	0.3	0.5	1.0	2.9
9:05 AM	0.1	0.5	0.2	0.2	1.2	3.5
9:33 AM	0.3	0.4	0.0	0.0	1.0	3.8
10:02 AM	0.5	0.1	0.5	0.3	1.0	3.7
10:31 AM	0.6	2.0	1.0	0.5	1.5	2.9
10:58 AM	3.0	3.1	2.3	2.2	2.0	0.0
11:35 AM	3.5	3.4	2.0	2.0	1.8	0.0
11:50 PM	Drug Dose A					

References

- Hoff JI, van Hilten BJ, and Roos RA (1999) A review of the assessment of dyskinesias. *Mov Disord* **14**, 737-43.
- Amanzio M, Monteverdi S, Giordano A, Soliveri P, Filippi P, and Geminiani G (2010) Impaired awareness of movement disorders in Parkinson's disease. *Brain Cogn* **72**, 337-46.
- Sitek EJ, Soltan W, Wieczorek D, Robowski P, Schinwelski M, and Slawek J (2011) Assessing self-awareness of dyskinesias in Parkinson's disease through movie materials. *Funct Neurol* **26**, 121-6.
- Stone AA, Shiffman S, Schwartz JE, Broderick JE, and Hufford MR (2003) Patient compliance with paper and electronic diaries. *Control Clin Trials* **24**, 182-99.
- Pechevis M, Clarke CE, Vieregge P, Khoshnood B, Deschaseaux-Voinet C, Berdeaux G, Ziegler M, and Trial Study G (2005) Effects of dyskinesias in Parkinson's disease on quality of life and health-related costs: a prospective European study. *Eur J Neurol* **12**, 956-63.
- Van Gerpen JA, Kumar N, Bower JH, Weigand S, and Ahlskog JE (2006) Levodopa-associated dyskinesia risk among Parkinson disease patients in Olmsted County, Minnesota, 1976-1990. *Arch Neurol* **63**, 205-9.
- Grandas F, Galiano ML, and Tabernerero C (1999) Risk factors for levodopa-induced dyskinesias in Parkinson's disease. *J Neurol* **246**, 1127-33.
- Carta M, Carlsson T, Munoz A, Kirik D, and Bjorklund A (2008) Serotonin-dopamine interaction in the induction and maintenance of L-DOPA-induced dyskinesias. *Prog Brain Res* **172**, 465-78.

Questions?

Contact Information:

Thomas Mera, MS
Sr. Biomedical Research Engineer
tmera@glneurotech.com
Great Lakes NeuroTechnologies

www.glneurotech.com/movement_disorders