

# Quantitative Parkinson's Gait Assessment: A high resolution measure of change in impairment

---

Tuesday July 22<sup>th</sup>, 2014

Starts at 12:00 PM EST

Presented by

Elizabeth Brokaw, PhD

## Outline

- Impairment due to Parkinson's Disease
- Deep Brain Stimulation
- Evaluations with wearable sensors
- DBS cessation research
- Role of Kinesia to improve gait outcomes

# Parkinson's Disease Motor Impairments


# Parkinson's Disease

- Wide range of motor symptoms
- Treatments outcomes often focus on the upper extremity
- Mobility is important for quality of life
- Increased fall risk after PD
  - 68.3% fell during one year



# Parkinson's Disease


- Gait
  - Freezing of gait
  - Slowed movement
  - Shortened stride length
  - Flat foot strike (shuffled steps)
  - Impaired balance and posture



# Deep Brain Stimulation

# Deep Brain Stimulation

- Became a standard treatment for PD in the 1990's
- Effective method for improving symptoms and reducing medication burden
- Typically indicated for
  - Tremor
  - Bradykinesia
  - Rigidity



# Deep Brain Stimulation for Gait Impairment

- Researcher have observed STN stimulation induced improvement in
  - Stride length
  - Walking speed
  - Freezing of gait
- Evaluation of settings and location is ongoing
  - Frequency effects
  - Pedunculopontine nucleus (PPN)




# Deep Brain Stimulation Response Time

## Deep Brain Stimulation

- Effects some symptoms quickly

Time For Effect on Motor Function After DBS Change  
(In Minutes)



- Effect on gait is slow and less predictable
  - Unknown final effect on gait and balance
  - Not optimized to improve gait and balance

# Wearable Sensors

# Wearable Sensors

- Quantitative Evaluation of Movement
  - Acceleration
  - Angular Velocity



## Benefits of Wearable Movement Sensors

- Objective measure of impairment
- High resolution
- Not confined to in clinic evaluations

# Quantifying Effect of DBS Cessation With the Kinesia Sensors


The goal is to examine changes in impairment related to changes in to DBS settings

# Kinesia

- Quantitative assessment of
  - Tremor
  - Bradykinesia
  - Dyskinesia
  - Gait
  - Freezing of gait
- For more information
  - <http://glneurotech.com/publications/>



# Kinesia Lower Extremity and Gait Evaluation



Published:

Heldman, D., Filipkowski, D. E., Riley, D. E., Whitney, C. M., Walter, B. L., Gunzler, S. a, Giuffrida, J.P. & Mera, T (2012). Automated motion sensor quantification of gait and lower extremity bradykinesia. *International conference of the IEEE EMBS. 2012.*

# Kinesia Sensor Placement

- Sensors placed on the more affected thigh, back and top of feet.





# Protocol

- Individuals with Parkinson's Disease and DBS
  - Started off medication and with DBS on
- Kinesia Evaluation: Unified Parkinson's Disease Rating Scale tasks

Toe Tapping



Leg Lifts



Arise



Gait



# Protocol

- Clinician UPDRS at study start DBS on and 3 hours after DBS off
- 2 Kinesia evaluations at each time point

Time After DBS Was Turned Off In Minutes



- Times of Kinesia evaluations
- Times of clinician evaluations

## Study Goals

- Evaluate changes in impairment over time
- Evaluate ability to minimize sensor number to reduce user burden
  - Sensors data from both legs
  - Sensor data from just the subject reported more affected limb

## Results

- 8 Individuals with Parkinson's Disease
- STN DBS implanted
- DBS surgery average of  $1.8 \pm 2.3$  years prior
- Average of  $14 \pm 1.5$  hours off medication

# Overall Effect of DBS Cessation

		Average Kinesia Score	Average Clinician Score
<b>Gait</b>	DBS on	$0.77 \pm 0.38$	$0.5 \pm 0.53$
	DBS off 3 hr	$1.00 \pm 0.45$	$0.88 \pm 0.99$
P value		<b>0.001 *</b>	0.28
<b>Toe Taps</b>			
	DBS on	$2.19 \pm 0.57$	$1.63 \pm 0.92$
	DBS off 3 hr	$2.58 \pm 0.49$	$2.38 \pm 1.06$
P value		<b>&lt;0.0001*</b>	<b>0.02*</b>
<b>Leg Lifts</b>			
	DBS on	$1.67 \pm 0.76$	$0.57 \pm 0.53$
	DBS off 3 hr	$2.54 \pm 0.94$	$2 \pm 0.58$
P value		<b>&lt;0.0001*</b>	<b>0.003*</b>

# Effect of DBS Cessation Over Time


## Gait

### Both Legs

#### Paired t-tests

- DBS On to 120 min  
( $p=0.002$ )

- 120 to 180 min  
( $p=0.53$ )



# Effect of DBS Cessation Over Time


## Gait

### More Impaired Leg

Paired t-tests

- DBS On to 120 min  
( $p=0.04$ )

- 120 to 180 min  
( $p=0.25$ )




# Effect of DBS Cessation Over Time

## Toe Tapping

### Both Legs

#### Paired t-tests

- DBS On to 15 min  
( $p=0.016$ )
- 15 to 180 min  
( $p=0.1$ )





# Effect of DBS Cessation Over Time


## Toe Tapping

### More Impaired Leg

Paired t-tests

- DBS On to 15 min  
( $p=0.02$ )

- 15 to 180 min  
( $p=0.13$ )




# Effect of DBS Cessation Over Time

## Leg Lifts

### Both Legs

#### Paired t-test

- DBS On to Off  
( $p=0.002$ )
- Off to 120 min  
( $p<0.001$ )
- 120 to 180 min  
( $p=0.4$ )




# Effect of DBS Cessation Over Time

## Leg Lifts

### More Impaired Leg


#### Paired t-test

- DBS On to Off ( $p=0.007$ )
- Off to 60 min ( $p<0.001$ )
- 60 to 180 min ( $p=0.99$ )



## Summary of DBS Cessation

- Increase in impairment after off 3 hour with DBS off
  - Except clinician gait score
- Different movements showed very different time response.



## Limitations

- Small sample size (study is ongoing)
- DBS cessation as a model for change in DBS settings


## Role of the Kinesia Systems



# Kinesia



Tune DBS settings  
in the clinic



Independent home  
assessments



# Kinesia for DBS


- Integrate remote evaluation and DBS tuning
  - This will improve knowledge of DBS effects
  - Allow for tuning of gait parameters





# Conclusions

- DBS changes over time



- Kinesia system
  - High resolution quantitative evaluation
  - Not limited to use in the clinic
  - Integration of Kinesia and DBS tuning could improve gait outcomes

# Acknowledgements



- Thomas Mera
- Dustin Heldman
- Joseph Giuffrida



- David Riley
- Benjamin Walter
- Steven Gunzler



- Alberto Espay
- Fredy Revilla

# Questions

For more information contact:

Elizabeth Brokaw

[Ebrokaw@glneurotech.com](mailto:Ebrokaw@glneurotech.com)