

Kinesia One Symptom Severity Report

Patient Name: Sample Patient
Patient ID: 103-129
Institution: Sample Hospital
Birthdate: 09/18/1949

Time	Rest Tremor	Postural Tremor	Kinetic Tremor	Finger Taps Speed	Finger Taps Amplitude	Finger Taps Rhythm	Hand Movements Speed	Hand Movements Amplitude	Hand Movements Rhythm	Rapid Alternating Movements Speed	Rapid Alternating Movements Amplitude	Rapid Alternating Movements Rhythm	Dyskinesia	
01/15/2015	9:52 AM	3.0	1.2	1.2	2.3	2.5	1.1	2.4	2.3	0.9	2.1	1.4	1.3	0.0
	12:52 PM	1.2	0.2	1.1	2.4	2.7	2.0	2.3	2.2	1.7	2.1	0.7	1.3	0.4
	Mean	2.1	0.7	1.1	2.4	2.6	1.5	2.4	2.3	1.3	2.1	1.1	1.3	0.2
	Fluctuation	0.9	0.5	0.1	0.0	0.1	0.5	0.1	0.1	0.4	0.0	0.3	0.0	0.2
01/31/2015	7:00 AM	3.1	2.3	1.6	2.0	1.6	0.9	1.9	1.0	0.8	1.9	0.9	0.9	0.0
	8:00 AM	3.4	2.8	1.5	1.9	1.6	1.2	1.7	0.9	1.0	1.7	0.9	0.9	0.0
	9:00 AM	3.4	1.4	1.4	1.7	1.2	1.2	1.7	0.8	0.6	1.7	0.9	1.1	0.0
	10:01 AM	3.4	2.2	1.7	1.9	1.6	1.4	1.7	1.0	0.9	1.7	0.6	0.7	0.0
	11:01 AM	3.2	1.6	1.6	1.5	0.3	0.6	1.6	0.5	0.7	1.8	0.7	0.8	0.0
	12:00 PM	0.0	0.0	1.2	1.6	0.9	1.1	1.4	0.6	0.7	1.7	0.6	0.4	1.4
	1:00 PM	2.6	0.6	1.3	1.8	1.5	1.0	1.7	1.3	1.0	1.6	0.5	0.6	0.0
	2:00 PM	2.3	0.7	1.6	2.0	2.0	1.1	1.7	1.4	1.0	1.8	0.7	0.8	0.0
	3:01 PM	2.9	1.6	1.4	1.7	2.0	0.8	1.6	1.8	0.8	1.7	0.6	0.6	0.0
	4:00 PM	3.5	1.9	1.5	1.8	1.4	1.2	1.8	1.4	1.0	1.8	0.8	0.7	0.0
	5:00 PM	3.1	1.5	1.6	1.6	1.2	1.2	1.5	1.2	1.0	1.7	0.8	0.7	0.0
	6:00 PM	3.3	1.6	1.8	1.7	1.3	1.4	1.7	1.2	1.2	1.8	0.8	0.7	0.0
	7:00 PM	2.9	1.4	1.8	1.8	1.7	1.5	1.6	1.1	0.9	1.9	0.8	1.1	0.0
	8:00 PM	3.0	1.0	1.5	1.9	1.2	1.3	1.7	1.3	1.2	1.9	0.8	0.7	0.0
	9:00 PM	3.0	1.5	1.7	1.7	0.8	1.4	1.4	0.4	0.5	1.8	0.7	0.8	0.0
	10:00 PM	2.9	0.8	1.3	2.0	1.8	1.0	1.6	0.6	0.6	1.9	0.7	0.6	0.0
	11:00 PM	3.2	0.0	1.5	1.6	0.8	0.9	1.8	1.1	1.0	1.9	0.7	1.0	0.7
Mean	2.9	1.3	1.5	1.8	1.3	1.1	1.6	1.0	0.9	1.8	0.7	0.8	0.1	
Fluctuation	0.8	0.7	0.2	0.1	0.5	0.2	0.1	0.4	0.2	0.1	0.1	0.2	0.4	

	Time	Rest Tremor	Postural Tremor	Kinetic Tremor	Finger Taps Speed	Finger Taps Amplitude	Finger Taps Rhythm	Hand Movements Speed	Hand Movements Amplitude	Hand Movements Rhythm	Rapid Alternating Movements Speed	Rapid Alternating Movements Amplitude	Rapid Alternating Movements Rhythm	Dyskinesia
02/10/2015	7:02 AM	0.0	0.0	1.0	1.9	1.7	1.3	1.7	1.4	0.9	1.3	0.2	0.6	1.1
	8:03 AM	0.0	0.0	1.0	1.9	1.8	1.1	1.7	1.9	1.0	1.5	0.4	0.6	1.2
	9:00 AM	1.3	0.0	1.1	1.8	2.2	0.9	1.7	1.7	1.0	1.4	0.3	0.6	0.2
	10:00 AM	0.0	0.0	1.1	1.8	2.3	1.0	2.0	2.3	1.0	1.4	0.5	0.6	1.2
	11:00 AM	2.3	0.0	1.2	1.6	1.6	0.9	1.6	1.3	1.2	1.4	0.5	0.7	0.2
	12:00 PM	2.7	0.0	1.1	1.8	1.9	0.8	1.5	1.4	0.8	1.6	0.7	0.6	0.4
	1:00 PM	2.5	0.0	1.2	1.8	1.9	0.9	1.7	1.7	0.8	1.5	0.6	0.4	0.2
	2:01 PM	0.0	0.0	1.1	1.7	1.7	0.8	1.7	1.6	1.7	1.4	0.5	0.5	1.0
	3:00 PM	0.0	0.0	1.5	1.8	1.7	0.7	1.5	1.5	1.5	1.3	0.4	0.5	1.8
	4:00 PM	0.0	0.0	1.1	1.6	1.9	0.9	1.7	1.9	1.1	1.5	0.3	0.5	1.9
	5:00 PM	0.0	0.0	1.2	1.5	1.5	0.8	1.6	1.8	0.9	1.5	0.6	0.2	1.7
	6:00 PM	2.3	0.1	1.3	1.5	1.8	1.2	1.7	1.7	1.1	1.5	0.6	0.8	0.4
	7:08 PM	2.4	0.2	1.0	1.6	1.9	0.8	1.8	1.8	1.1	1.5	0.7	0.5	0.3
	8:00 PM	0.0	0.0	1.0	1.6	2.1	1.5	1.3	1.1	0.7	1.3	0.6	0.7	1.5
	9:00 PM	0.0	0.0	1.5	1.8	2.3	0.7	1.5	1.5	1.0	1.2	0.4	0.3	1.8
11:01 PM	0.0	0.0	1.4	1.8	2.0	0.7	1.9	1.7	1.0	1.5	0.7	0.6	1.6	
	Mean	0.8	0.0	1.2	1.7	1.9	0.9	1.7	1.6	1.0	1.4	0.5	0.5	1.0
	Fluctuation	1.1	0.0	0.2	0.1	0.2	0.2	0.2	0.3	0.3	0.1	0.1	0.2	0.6
02/25/2015	7:02 AM	0.0	0.1	1.0	1.8	1.9	0.9	1.6	1.9	0.7	1.5	0.7	0.5	1.4
	8:01 AM	0.0	0.0	1.2	1.7	2.1	1.0	1.8	2.3	0.9	1.1	0.7	0.5	1.2
	9:00 AM	0.0	0.0	1.2	1.9	2.1	0.8	1.7	2.0	0.8	1.5	0.5	0.9	1.4
	10:01 AM	0.0	0.0	1.0	2.0	2.5	1.2	1.7	2.0	1.1	1.3	0.5	0.6	1.0
	11:00 AM	2.5	0.6	1.2	1.9	1.9	1.6	1.8	2.0	1.0	1.6	0.9	0.9	0.0
	12:00 PM	2.0	0.0	1.3	1.8	2.0	0.9	1.8	2.0	0.7	1.7	0.9	0.6	0.4
	1:00 PM	2.1	0.0	1.3	1.8	2.0	1.1	1.7	1.9	1.0	1.6	0.8	1.0	0.6
	2:01 PM	3.0	1.1	1.4	1.6	1.9	1.0	2.0	2.3	1.2	1.7	1.2	0.5	0.0
	3:01 PM	2.2	0.8	1.3	1.9	2.0	1.3	1.9	2.2	0.9	1.6	0.9	0.6	0.0
	4:00 PM	2.1	0.6	1.2	2.1	2.3	0.9	1.6	2.0	1.0	1.6	1.1	0.6	0.0
	5:00 PM	1.6	0.3	1.3	1.9	2.4	1.3	1.7	2.0	0.6	1.6	0.8	0.6	0.3
	6:00 PM	1.4	0.2	1.7	1.5	2.0	1.0	1.7	2.1	0.8	1.5	0.8	0.6	0.3
	7:01 PM	0.0	0.0	1.0	1.7	2.1	1.3	1.5	1.6	0.8	1.4	0.6	0.3	1.2
	8:00 PM	0.0	0.0	1.3	1.6	2.1	1.0	1.8	2.0	1.1	1.5	0.8	0.4	1.5
	9:00 PM	0.0	0.0	1.2	1.8	2.2	0.7	1.6	1.7	0.8	1.5	0.6	0.6	1.2
10:00 PM	0.0	0.0	1.1	1.9	2.3	1.1	1.7	1.9	0.9	1.5	0.9	0.5	1.6	
11:01 PM	0.9	0.0	1.3	1.8	2.1	0.8	1.8	2.1	0.8	1.6	0.7	0.5	0.4	
	Mean	1.0	0.2	1.2	1.8	2.1	1.0	1.7	2.0	0.9	1.5	0.8	0.6	0.7
	Fluctuation	1.1	0.3	0.2	0.1	0.2	0.2	0.1	0.2	0.2	0.1	0.2	0.2	0.6